

MET'S Institute of Pharmacy

ACTIVITY REPORT

Name of Activity - Installation, demonstration and training to student on Lyophilizer / Freeze Drier LAIPL-LI-10-MR-80

Day and Date- Thursday 05/02/2026

	Staff	Student	Class and Course
Participants	03	16	M. Pharm (QAT)

Details of the Program

A technical representative from Labindia Analytical Instruments Pvt. Ltd. visited the Institute of Pharmacy (IOP) for the installation, demonstration, and operational training of the Lyophilizer/Freeze Dryer (Model: LAIPL-LI-10-MR-80). The session included equipment setup, working principles, operating procedures, safety precautions, routine maintenance, and troubleshooting guidelines. Students and staff of the QAT department actively participated in the hands-on training and gained practical exposure to the instrument's functionality and applications.

Programme outcomes attained

PO1	Pharmacy knowledge	✓	PO7	Pharmaceutical ethics	✓
PO2	Planning abilities	✓	PO8	Communication	✓
PO3	Problem analysis		PO9	The pharmacist and society	✓
PO4	Modern tool usage	✓	PO10	Environment and sustainability	✓
PO5	Leadership skills	✓	PO11	Lifelong learning	✓
PO6	Professional identity	✓			

Objective of activity	1. Practical experience with Factory Acceptance Testing 2. Understanding quality assurance processes in pharmaceutical equipment manufacturing
Benefits of the activity	The activity provided hands-on experience, industry exposure, and practical learning for students, enhancing their understanding and employability in quality assurance.
Social Media links	

MET'S Institute of Pharmacy

Photos of the event:




Dr. S. J. Kshirsagar

Principal

Dr. R. S. Ahire 
Dr. E. D. Ahire 
Ms. P. M. Gadakh 

coordinators


Dr. G. S. Deokar

Coordinator & HOD
(GA)